



WARRANTY INSPECTION AND MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups according to the warranty inspection schedule below.

High patient weight and/or impact level may require more frequent inspections. We recommend you visually inspect the following applicable parts for excessive wear and fatigue at each warranty inspection.

- Composites and Adapters
- CPI Sock
- Foot Shell

WARRANTY INSPECTION SCHEDULE FOR CELSUS: ANNUALLY.

TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 AM – 5:30 PM (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.



**COLLEGE PARK
INDUSTRIES, INC.**
27955 College Park Dr.
Warren, MI 48088 USA

©2018 College Park Industries, Inc. All rights reserved. College Park, Celsus, Intelliweave, and Gait Matching are registered trademarks of College Park Industries, Inc.


EMERGO EUROPE
Prinsessegracht 20, 2514 AP
The Hague, Netherlands

CE
MADE IN THE USA
126 INS CS TIS 171206

**Australian Sponsor
EMERGO AUSTRALIA**
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia



*technical
instructions*



This diagram is to help familiarize you with the unique parts of the Celsus. These parts are referenced in the instructions and used when speaking with a technical service representative.

KEY COMPONENTS

1. Integrated Pyramid
2. Housing (Ankle)
3. Fulcrum Pad
4. Toe Spring
5. Heel Wedge (permanent)
6. Heel Spring

PACKAGE CONTENTS

- (1) Celsus foot
- (1) Foot Shell
- (1) CPI sock

TOOLS RECOMMENDED

- 4mm Hex Key
- FootHorn

⚠ WARNING

- Do not expose this product to corrosive materials, salt water or pH extremes.
- Contaminants such as dirt and the use of lubricants or powder, may effect the function of the CPI Sock and lead to noise.
- Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.

⚠ PROTECTIVE COVER ON DOME

Remove the protective cover on dome after alignment is completed and before patient leaves clinic.

MOUNTING

Use only high quality endoskeletal components.

ASSEMBLY AND DISASSEMBLY (for sock replacement)

Use the FootHorn to don and doff the foot shell. Remove the CPI Sock and replace as needed. Any further disassembly or modification of components will void the warranty.

STATIC ALIGNMENT

For optimal function have the patient's weight balanced between the forefoot and heel. The load line divides the foot at the 1/3 heel lever and 2/3 toe lever.

DYNAMIC ADJUSTMENTS

DESIRED RESULT	ALIGNMENT CHANGE
Firmer Toe Response	Plantarflex the Celsus or move load line posterior
Softer Toe Response	Dorsiflex the Celsus or move load line anterior
Firmer Heel Response	Dorsiflex the Celsus or move load line anterior
Softer Heel Response	Plantarflex the Celsus or move load line posterior

GAIT MATCHING® GUIDELINES

Contact College Park Technical Service if you have questions about modifying the gait match.