

## WARRANTY INSPECTION AND MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups according to the warranty inspection schedule below.

High patient weight and/or impact level may require more frequent inspections. We recommend you visually inspect the following applicable parts for excessive wear and fatigue at each warranty inspection.

- Composites and Adapters
- CPI Sock
- Foot Shell

**WARRANTY INSPECTION SCHEDULE FOR VELOCITY: ANNUALLY.**

## TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 AM – 5:30 PM (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.



**COLLEGE PARK  
INDUSTRIES, INC**  
27955 College Park Dr.  
Warren, MI 48088 USA

**EC REP EMERGO EUROPE**  
Prinsessegracht 20, 2514 AP  
The Hague, Netherlands


**CE**  
**MADE IN THE USA**  
110 INS VLT IS 171211

**Australian Sponsor  
EMERGO AUSTRALIA**  
Level 20, Tower II  
Darling Park  
201 Sussex Street  
Sydney, NSW 2000  
Australia

©2018 College Park Industries, Inc. All rights reserved. College Park, Velocity, Gait Matching, and Integrated Spring Technology (IST) are registered trademarks of College Park Industries, Inc. US Patent No. D655009.



*technical  
instructions*



This diagram is to help familiarize you with the unique parts of the Velocity. These parts are referenced in the instructions and used when speaking with a technical service representative.

## KEY COMPONENTS

1. Integrated Pyramid
2. Housing
3. Posterior Toe Plate
4. Anterior Toe Plate
5. Posterior Toe
6. Anterior Toe
7. Foreheel
8. Heel Pad
9. Toe Pad
10. Heel Wedge (optional)
11. Toe Wedge (optional)

## PACKAGE CONTENTS

- (1) Velocity Foot
- (1) Foot Shell
- (1) CPI Sock
- (1) Velocity Wedge Kit

## TOOLS REQUIRED

- (1) 4mm Hex Key

## TOOLS RECOMMENDED

- Foot Horn

## ⚠️ PROTECTIVE COVER ON DOME

Remove the protective cover on dome after alignment is completed and before patient leaves clinic.

## MOUNTING

Use only high quality proximal endoskeletal components.

## ASSEMBLY AND DISASSEMBLY (for sock replacement)

Use the FootHorn to don and doff the foot shell. Remove the CPI Sock and replace as needed. Any further disassembly or modification of components will void the warranty.

## STATIC ALIGNMENT

For optimal function have the patient's weight balanced between the forefoot and heel. The load line divides the foot at 1/3 heel lever.

## DYNAMIC ADJUSTMENTS

DESIRED RESULT	ALIGNMENT CHANGE	COMPONENT CHANGE
<b>Firmer Toe Response</b>	Plantarflex the Velocity or move load line posterior	Install Toe Wedge*
<b>Softer Toe Response</b>	Dorsiflex the Velocity or move load line anterior	Remove Toe Wedge*
<b>Firmer Heel Response</b>	Dorsiflex the Velocity or move load line anterior	Install Heel Wedge*
<b>Softer Heel Response</b>	Plantarflex the Velocity or move load line posterior	Remove Heel Wedge*

\* To install or remove toe and heel wedges, refer to the Velocity Wedge Kit Instructions. Pre-installed toe wedges cannot be removed.

## GAIT MATCHING® GUIDELINES

Contact College Park Technical Service if you have questions about modifying the gait match.

## ⚠️ WARNING

- Do not expose this product to corrosive materials, salt water or pH extremes.
- Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.