

Regal Silicone Prostheses

Regal silicone prosthetic gloves - Size measurement chart

For Model 101L, 102LS, 102L, 101LL, 102LSL, 102LL (wrist or above wrist amputation)

Male - gloves thickness: 2-2.5mm										Female - gloves thickness: 2-2.5mm																			
OttoBock gloves	Regal gloves		(mm)				Standard length	Long length	Catalogue 2015 page		OttoBock gloves	Regal gloves		(mm)				Standard length	Long length	Catalogue 2015 page									
			C1	C3	C11	L5								E	F	C1	C3				C11	L5	E	F					
8S4=238x92L	XL	Left	235	196	80	74	310	544	P.54	Adult	8S5=208x89L	XL	Left	220	183	79	76	314	535	P.64									
8S4=244x94R	XL	Right	239	210	84	76	310	538	P.54		8S5=210x89R	XL	Right	223	183	76	80	323	536	P.64									
8S4=228x88L	XL1	Left	233	196	82	73	395	538	P.56		-	L	Left	205	179	72	67	297	487	P.82									
8S4=228x84L		8S4=221x81L	-	L	Right	230	194	77	77		373	530	P.56	-	L	Right	206	179	71	65	290	487	P.82						
8S4=228x89R	XL1	Right	230	194	77	77	373	530	P.56		8S5=195x78L	L1	Left	202	170	66	74	315	532	P.70									
8S4=225x82R		8S4=222x84R	8S5=195x78R	L1	Right	203	174	64	71		319	532	P.70																
8S4=232x94L	L	Left	230	190	74	78	320	533	P.58		8S5=195x79L	L2	Left	202	174	72	70	313	470	P.84									
8S4=230x93R	L	Right	231	196	75	77	318	530	P.58		8S5=200x79R	L2	Right	203	174	72	67	310	460	P.84									
-	L1	Left	220	193	76	67	335	556	P.60																				
-	L1	Right	233	200	81	66	336	550	P.60																				
-	L2	Left	228	187	73	76	370	545	P.62																				
-	L2	Right	229	197	77	76	373	546	P.62																				
8S4=220x91L	LM	Left	220	183	79	76	314	535	P.64		8S5=194x82L	LM	Left	189	163	62	67	288	481	P.76									
8S4=218x85L											8S5=192x78L																		
8S4=214x82L																													
8S4=213x85L											8S4=212x93L																		
8S4=218x85R	LM	Right	223	183	76	80	323	536	P.64		8S5=191x78R	LM	Right	198	166	64	70	295	488	P.76									
8S4=218x83R											8S5=189x84R																		
8S4=215x83R																													
8S4=214x90R											8S4=215x93R																		
8S4=211x88L	M	Left	209	185	72	70	327	542	P.66		-	LM1	Left	190	162	63	70	290	501	P.78									
8S4=208x85L											8S4=207x86L																		
8S4=207x86L											8S4=206x87L																		
8S4=212x86R											8S4=212x83R																		
8S4=212x83R	M	Right	212	188	72	72	320	532	P.66		-	LM1	Right	198	162	65	71	303	497	P.78									
8S4=215x86R											8S4=209x86R																		
8S4=206x80L											8S4=205x81L																		
8S4=205x81R											8S4=205x80R																		
8S4=203x85L	M2	Left	202	170	66	74	315	532	P.70		8S5=182x84L	M	Left	182	156	66	69	293	452	P.86									
8S4=203x83L											8S5=180x80L																		
8S4=206x85R										8S4=199x82R																			
8S4=202x74L										8S4=206x76R																			
8S4=202x74L	M3	Left	213	176	67	73	330	528	P.72	8S5=190x93L	M	Right	184	155	66	68	293	450	P.86										
8S4=206x85R										8S5=188x79L																			
8S4=199x82R										8S5=190x93R	MS	Left	188	162	62	69	307	463	P.80										
8S4=202x74L										8S5=190x77R	MS	Right	190	160	64	70	317	463	P.80										
8S4=206x76R	8S5=188x79R	-	MS2	Left	181	157	61	65	335	526										P.90									
-	M4	Left	210	177	71	69	344	512	P.74	-	MS2	Right	185	157	63	64	330	522	P.90										
-										8S5=185x75L	S	Left	180	160	62	60	286	479	P.88										
-										8S5=184x78L																			
-										8S5=184x75L																			
-	M4	Right	202	167	67	67	340	530	P.74	8S5=187x74R	S	Right	188	162	62	62	292	468	P.88										
-										8S5=184x78R																			
-										8S5=181x75R																			
-										MS	Left	189	163	62	67	288	481	P.76	8S5=174x74L	S1	Left	178	144	54	65	299	481	P.92	
-	8S5=165x72L																												
-	MS	Right	198	166	64	70	295	488	P.76										8S5=175x76R	S1	Right	175	147	54	65	298	477	P.92	
-																			8S5=167x72R										
-	MS1	Left	190	162	63	70	290	501	P.78	8S6=170x65L	CL	Left	178	144	54	65	299	481	P.92										
-										8S6=170x65R										CL	Right	175	147	54	65	298	477	P.92	
-										MS1	Right	198	162	65	71	303	497	P.78											
-																			S	Left	188	162	62	69	307	463	P.80		
-	S	Right	190	160	64	70	317	463	P.80																				
-										CL	Left	188	162	62	69	307	463	P.80											
-	CL	Left	188	162	62	69	307	463	P.80																				
-										CL	Right	190	160	64	70	317	463	P.80											
-										CL1	Left	181	157	61	65	335	526	P.90											
-																			CL1	Right	185	157	63	64	330	522	P.90		
8S6=170x65L	CM	Left	178	144	54	65	299	481	P.92	8S6=168x70L	CM	Left	173	142	57	62	310	473	P.94										
8S6=168x70L										8S6=165x68L																			
8S6=165x68L																													
8S6=170x65R										CM	Right	175	147	54	65	298	477	P.92	8S6=166x70R	CM	Right	165	138	55	64	295	473	P.94	
8S6=166x70R	CS	Left	141	119	51	48	255	443	P.96	8S6=158x68R	CS	Left	141	119	51	48	255	443	P.96										
8S6=158x68R										8S6=142x50L																			
8S6=142x50L										8S6=138x58L																			
8S6=138x58L										8S6=139x51R	CS	Right	140	123	50	49	245	437	P.96										
8S6=139x51R	8S6=138x58R																												
-	CX	Left	125	115	43	45	214	391	P.98	-	CX	Left	125	115	43	45	214	391	P.98										
-										-										CX	Right	131	120	43	44	223	383	P.98	

The above information is for reference only. Please contact our technicians if you are not sure which size to choose.