


AK & BK
KITS



BELOW KNEE KITS

BK KIT #1

CONTENTS:

- Cushion liner
- Suspension sleeve
- Expulsion valve
- 4 ear adaptor
- Short tube
- Tube clamp
- College Park Breeze foot
- BK foam block
- BK cosmetic stockings



BK KIT #2

CONTENTS:

- Cushion liner
- Suspension sleeve
- Expulsion valve
- 4 ear adaptor
- Short tube
- Tube clamp
- College Park Celsus foot
- BK foam block
- BK cosmetic stockings



BK KIT #3

CONTENTS:

- 4 ear adaptor
- Short tube
- Tube clamp
- SACH adaptor
- SACH foot
- BK foam block
- BK cosmetic stockings



BK KIT #4

CONTENTS:

- Liner with pin attachment
- Shuttle lock
- Short tube
- Tube clamp
- SACH adaptor
- SACH foot
- BK foam block
- BK cosmetic stockings



BK KIT #5

CONTENTS:

- 4 ear adaptor
- Short tube
- Tube clamp
- Single Axis adaptor
- Single Axis foot
- BK foam block
- BK cosmetic stockings



BK KIT #6

CONTENTS:

- Liner with pin attachment
- Shuttle lock
- Short tube
- Tube clamp
- Single Axis adaptor
- Single axis foot
- BK foam block
- BK cosmetic stockings



BK KIT #7

CONTENTS:

- Liner with pin attachment
- Shuttle lock
- Short tube
- Tube clamp
- College Park Breeze foot
- BK foam block
- BK cosmetic stockings




BK KIT #8

CONTENTS:

- Liner with pin attachment
- Shuttle lock
- Short tube
- Tube clamp
- College Park Celsus foot
- BK foam block
- BK cosmetic stockings





ABOVE KNEE KITS

AK KIT #1



CONTENTS:

- 3 prong rotating receiver
- Safety Knee
- Long tube
- Tube clamp
- SACH adaptor
- SACH foot
- AK foam cover
- AK cosmetic stockings

AK KIT #2



CONTENTS:

- 3 prong rotating receiver
- Safety Knee
- Long tube
- Tube clamp
- Single Axis adaptor
- Single Axis foot
- AK foam cover
- AK cosmetic stockings

AK KIT #3



CONTENTS:

- Liner with pin attachment
- Shuttle lock
- A305 knee joint
- Long tube
- College Park Breeze foot
- AK foam cover
- AK cosmetic stockings

AK KIT #4



CONTENTS:

- Liner with pin attachment
- Shuttle lock
- Ti Habermann type knee joint
- Tube clamp
- Long tube
- College Park Breeze foot
- AK foam cover
- AK cosmetic stockings

AK KIT #5



CONTENTS:

- Liner with pin attachment
- Shuttle lock
- 6H25 Pneumatic knee joint
- Long tube
- College Park Celsus foot
- AK foam cover
- AK cosmetic stockings

AK KIT #6



CONTENTS:

- Liner with pin attachment
- Shuttle lock
- 6H30 hydraulic knee joint
- Long tube
- College Park Celsus foot
- AK foam cover
- AK cosmetic stockings