

AK & BK KITS



# BELOW KNEE KITS

# BK Kit #1

| Silicone cushion liner | Suspension sleeve | College Park Breeze foot |

## CONTENTS:

- Wagner silicone cushion liner
- Wagner suspension sleeve
- Expulsion valve
- 4 ear adaptor
- Short tube
- Tube clamp
- College Park Breeze foot
- BK foam block
- Cosmetic stockings



# BK Kit #2

| Silicone cushion liner | Suspension sleeve | College Park Celsus foot |



## CONTENTS:


- Wagner silicone cushion liner
- Wagner suspension sleeve
- Expulsion valve
- 4 ear adaptor
- Short tube
- Tube clamp
- College Park Celsus foot
- BK foam block
- Cosmetic stockings

## BK Kit #3

| Basic kit | SACH foot |

### CONTENTS:

- 4 ear adaptor
- Short tube
- Tube clamp
- SACH foot adaptor
- SACH foot
- BK foam block
- Cosmetic stocking



## BK Kit #4

| Silicone liner | Shuttle lock | SACH foot

### CONTENTS:

- Wagner silicone liner
- Shuttle lock
- Short tube
- Tube clamp
- SACH foot adaptor
- SACH foot
- BK foam block
- Cosmetic stocking



# BK Kit #5

| Basic kit | Single Axis foot |

## CONTENTS:

- 4 ear adaptor
- Short tube
- Tube clamp
- Single axis foot adaptor
- Single axis foot
- BK foam block
- Cosmetic stockings



# BK Kit #6

| Silicone liner | Shuttle lock | Single Axis foot

## CONTENTS:

- Wagner silicone liner
- Shuttle lock
- Short tube
- Tube clamp
- Single axis foot adaptor
- Single axis foot
- BK foam block
- Cosmetic stockings



# BK Kit #7

| Silicone liner | Shuttle lock | College Park Breeze foot |

## CONTENTS:

- Wagner silicone liner
- Shuttle lock
- Short tube
- Tube clamp
- College Park Breeze foot
- BK foam block
- Cosmetic stockings




# BK Kit #8

| Silicone liner | Shuttle lock | College Park Celsus foot |

## CONTENTS:

- Wagner silicone liner
- Shuttle lock
- Short tube
- Tube clamp
- College Park Celsus foot
- BK foam block
- Cosmetic stocking





# ABOVE KNEE KITS

# AK Kit #1

| Basic kit | SACH foot |

## CONTENTS:

- 3 prong rotating receiver
- Safety knee
- Long tube
- Tube clamp
- SACH foot adaptor
- SACH foot
- AK foam cover
- Cosmetic stockings



# AK Kit #2

| Basic kit | Single axis foot |

## CONTENTS:

- 3 prong rotating receiver
- Safety knee
- Long tube
- Tube clamp
- Single axis foot adaptor
- Single axis foot
- AK foam cover
- Cosmetic stockings





## AK Kit #3

| Silicone liner | Shuttle lock | A305 knee | College Park Breeze foot |

### CONTENTS:

- Wagner silicone liner
- Shuttle lock
- A305 knee joint
- Long tube
- College Park Breeze foot
- AK foam cover
- Cosmetic stockings



## AK Kit #4

| Silicone liner | Shuttle lock | Habermann type knee | College Park Breeze foot |

### CONTENTS:

- Wagner silicone liner
- Shuttle lock
- Ti Habermann type knee joint
- Tube clamp
- Long tube
- College Park Breeze foot
- AK foam cover
- Cosmetic stockings




# AK Kit #5

| Silicone liner | Shuttle lock | Pneumatic Knee Joint | College Park Celsus foot |

## CONTENTS:

- Wagner silicone liner
- Shuttle lock
- 6H25 pneumatic knee joint
- Long tube
- College Park Celsus foot
- AK foam cover
- Cosmetic stockings





12 Dahlia Street  
Northmead  
Benoni, South Africa  
[www.xprosortho.com](http://www.xprosortho.com)

T: +27 11 849 7650  
E: [orders.xpros@gmail.com](mailto:orders.xpros@gmail.com)  
[office.xpros@gmail.com](mailto:office.xpros@gmail.com)