

ORTHOTICS CATALOGUE 2023



CONTENTS

XPROS ORTHOTICS CATALOGUE 2023

BACK BRACES

| 01 |

NECK COLLARS

| 04 |

ARM BRACES

| 06 |

KNEE BRACES

| 08 |

WALKING AIDS

| 11 |

COMPRESSION STOCKINGS


| 13 |

ORTHOPAEDIC BOOTS

| 15 |

MATERIALS

| 17 |



BACK BRACES



TLISO BRACE

Used to limit movement in the thoracic, lumbar and sacral regions of the spine. Absorptive padding for comfort, easy adjustment and universal sizing.

INDICATIONS:

- Arthritic back pain
- Spondylolysis
- Herniated disc
- Thoracolumbar muscle weakness
- Universal sizing



TLISO BRACE

Used to limit movement in the thoracic, lumbar and sacral regions of the spine. Padded lumbar support, easy adjustment pulley system. Y bar provides frontal height options .

INDICATIONS:

- Herniated disc and compression fractures
- Scoliosis
- Post-op spinal fusions
- Sizes: S - XXL



SOMI BRACE

Lightweight rigid cervical brace aligning the neck and spine. Immobilisation of cervical spine prevents rotation as well as A/P and M/L movement.

Fully adjustable setting - flexion, neutral or hyper-extension. Universal sizing with removable chin plate.

INDICATIONS:

- Compression fractures
- Severe ligament or muscle strains
- Post-op support
- Universal sizing



LUMBAR SACRAL ORTHOSIS

Low profile immobilising orthosis for the lumbar spine. Discreet fit with a contour hugging design. Easily adjustable straps for a comfortable fit.

INDICATIONS:

- Acute lumbar spine discomfort
- Support of lower spinal ligaments/muscles
- Post-op lumbar immobilisation
- Sizing: S - XXL



LUMBAR SACRAL CORSET

A soft lower back support brace made from breathable elastic material, with 4 metal stays and velcro straps.

INDICATIONS:

- Acute and chronic pain
- Post-op immobilisation
- Sizing: S - XXXL

S	M	L	XL	XXL	XXXL
75cm 85cm	85cm 95cm	95cm 105cm	105cm 115cm	115cm 125cm	125cm 135cm



NECK COLLARS



SOFT COLLAR

A soft cervical collar used to provide temporary neck support and restrict movement. Made from a lightweight, breathable material.

INDICATIONS:

- Acute neck pain
- Sizing: S - L



PHILADELPHIA COLLAR

A semi-rigid cervical collar used to prevent head and neck movement. Made from lightweight plastazote with tracheal opening.

INDICATIONS:

- Cervical hernia
- Muscle and ligament strains
- Osteoarthritis



RIGID CERVICAL COLLAR

A lightweight rigid cervical collar which alleviates pain and contraction by immobilising the cervical spine and restricting movement. Universal sizing, with height adjustable supports and tracheal opening.

INDICATIONS:

- Cervical hernia
- Cervical fracture
- Post-op immobilisation

ARM BRACES



THUMB SPIKA WRIST BRACE

Anatomically-contoured aluminium wrist stay with thumb immobilised in the neutral position. Wrist straps reduce movement.

INDICATIONS:

- Carpal tunnel syndrome
- Abductor/flexor pollicis strains
- Post-op immobilisation



ARM IMMOBILISER

Fabric arm cradle with over shoulder strap.

INDICATIONS:

- Provides immobilisation for fractures of the ulna and radius



ELBOW ROM BRACE

A post-operative elbow brace that restricts and limits movement of the elbow.

INDICATIONS:

- Post-op protection and immobilisation of elbow
- Velcro straps used to fasten the brace
- Adjustable degree of flexion and extension 10 to 90 degrees
- Available in left and right



KNEE BRACES



HINGED KNEE BRACE

Made from neoprene material and anatomically designed. Steel joints on both sides provide maximum support to allow controlled movement of the joint.

INDICATIONS:

- Ligament injuries and instability
- Meniscus lesions
- First or second-degree strain
- Post-op protection and immobilisation

S	M	L	XL	XXL	XXXL
31cm 34cm	34cm 37cm	37cm 40cm	40cm 44cm	44cm 48cm	48cm 52cm



CRISS CROSS KNEE BRACE

Made from neoprene material and anatomically designed. Padded patella support with cross straps, provides anterior cruciate ligament and patella stabilisation.

INDICATIONS:

- Cruciate ligament injuries and knee instability
- Patella stabilisation
- Meniscus lesions
- Provides post-operative protection and stabilisation

S	M	L	XL	XXL	XXXL
31cm 34cm	34cm 37cm	37cm 40cm	40cm 44cm	44cm 48cm	48cm 52cm



PATELLA SUPPORT KNEE BRACE

Made from neoprene material and anatomically designed. Open patella with support pad.

INDICATIONS:

- Arthritis, bursitis, tendinitis
- Strain
- Meniscus lesions
- Ligament instability
- Sports injuries
- Provides post-operative protection

S	M	L	XL	XXL	XXXL
31cm 34cm	34cm 37cm	37cm 40cm	40cm 44cm	44cm 48cm	48cm 52cm



KNEE ROM BRACE TELESCOPIC

A post-operative knee brace that restricts and limits movement of the knee joint. The telescopic stays can be adjusted to accommodate different heights, from 20" to 26".

INDICATIONS:

- Post-op protection and immobilisation of knee joint
- Velcro straps used to fasten the brace
- Adjustable degree of flexion and extension 10 to 90 degrees



KNEE ROM BRACE PAEDIATRIC

A post-operative knee brace for children that restricts or limits movement of the knee joint.

INDICATIONS:

- Post-op protection and immobilisation of knee joint
- Velcro straps used to fasten the brace
- Adjustable degree of flexion and extension 10 to 90 degrees



ANKLE SUPPORT WITH MALLEOLUS GEL PAD

Made from neoprene material. Malleolus supported by anatomically designed gel pads.

Supports the ankle, promotes pain relief and reduces edema.

INDICATIONS:

- Acute ligament injuries and instability
- Arthritis and osteoarthritis
- Post-op protection and support

S	M	L	XL	XXL	XXXL
21cm 23cm	23cm 25cm	25cm 27cm	27cm 29cm	29cm 31cm	31cm 33cm



WALKING AIDS



ELBOW CRUTCHES

Used to support patients with mobility issues and lower extremity injuries.

FEATURES:

- Adjustable height range (70.3cm - 92.7cm)
- Weight limit: 120kg



WALKING FRAME

Used by patients recovering from leg or back injuries as well as hip and knee replacements. Also used by patients with mobility and balance issues.

FEATURES:

- Foldable with 2 push button release system
- Adjustable height range (80cm - 87cm)
- Weight limit: 110kg



FOLDING ROLLATOR

Used by patients recovering from leg or back injuries as well as hip and knee replacements. The two front wheels enable energy efficient propulsion, while maintaining a stable base of support.

FEATURES:

- Foldable with 2 push button release system
- Adjustable height range (80cm - 87cm)
- Weight limit: 110kg



COMPRESSION STOCKINGS



KNEE HIGH STOCKINGS

140 DEN | CLASS 1 | CLASS 2

S - XXXL

FEATURES:

- Option of: open toe / closed toe
- Non-ladder fabric
- Comfortable fit
- Reinforced heel and toe band



THIGH HIGH STOCKINGS

140 DEN | CLASS 2

S - XXXL

FEATURES:

- Option of: open toe / closed toe
- Non-ladder fabric
- Comfortable fit
- Reinforced heel and toe band



PANTYHOSE

140 DEN | CLASS 2

S - XXXL

FEATURES:

- Option of: open toe / closed toe
- Non-ladder fabric
- Comfortable fit
- Reinforced heel and toe band



ORTHOPAEDIC BOOTS

AFO BOOTS

Quality boot that opens fully to easily accommodate an AFO.



FEATURES:

- Lightweight design with velcro straps
- Removable orthopaedic insole
- Wide toe box
- Rear and forward opening facilitates donning

SPECIFICATIONS:

Child sizes:	5 - 12
Youth sizes:	1 - 5
Adult sizes:	6 - 7

LTT BOOTS

FEATURES:

- Reinforced back strap
- Heel and toe counters
- High ankle lacing for stability and support
- Maximum opening for easy foot insertion

SPECIFICATIONS:

Child sizes:	UK3 - UK13
Youth sizes:	UK1 - UK5
Adult sizes:	UK6 - UK12
Colour:	Black



DENIS BROWNE BAR & TP BOOTS

Lasted shoes used for the treatment of Talipes Equinovarus and used with a DB bar



FEATURES:

- Rigid shank and outer sole
- Posterior inspection port

SPECIFICATIONS:

Sizes:	0000 - 5
--------	----------





MATERIALS



EVA LINING

ORTHOTIC LINING MATERIAL

SPECIFICATIONS:

Sheet size:	1.5m x 1m
Thickness:	1.5mm
Colour:	Teal Grey Multi




DIABETIC INNERSOLE LINING

SEMI-PERFORATED MICROFIBER LINING
+
2MM PORON LAYER

SPECIFICATIONS:

Sheet size:	11.5cm x 6m
Thickness:	2mm poron
Colour:	Beige




PELITE

30 SHORE

SOFT SHORE LINING MATERIAL FOR ORTHOTICS

SPECIFICATIONS:

Sheet size:	2m x 1m
Thickness:	3mm 5mm 10mm
Colour:	White
Shore:	30




MEMOFLEX

40 - 45 SHORE

PADDING MATERIAL FOR PROSTHETIC SOCKETS
AND ORTHOTIC DEVICES

SPECIFICATIONS:

Sheet size:	1550mm x 720mm
Thickness:	5mm
Colour:	Taupe Brown Black
Shore:	40-45



MEMOFLEX

50 - 55 SHORE

PADDING MATERIAL FOR PROSTHETIC SOCKETS
AND ORTHOTIC DEVICES

SPECIFICATIONS:

Sheet size:	1425mm x 660mm
Thickness:	3mm 5mm 10mm
Colour:	Taupe Black
Shore:	50-55



SUPERIOR SOLE

DURABLE, LIGHTWEIGHT MATERIAL USED FOR SHOE RAISES

SPECIFICATIONS:

Sheet size:	1850mm x 1200mm
Thickness:	5mm 10mm 15mm
Colour:	Black White
Shore:	50-55



46 6th Street
Northmead
Benoni, South Africa

T: +27 11 849 7650

E: orders.xpros@gmail.com

www.xprosortho.com