


## Regal silicone prosthetic gloves - Size measurement chart

Measure Palm Circumference **C1** of the sound hand then choose a few closest Sizes from Regal Size Measurement Chart. From the Sizes chosen, select the appearance that closely resemble to the sound hand.

Note: some sizes are unisex.



How to measure : See page 139

### For Model 103 (partial hand amputation)

Male - gloves thickness: 2-2.5mm								Catalogue 2015 page	Female - gloves thickness: 2-2.5mm								Catalogue 2015 page
Regal gloves		(mm)					Regal gloves		(mm)								
		C1	C3	C11	L5	E			C1	C3	C11	L5	E				
XL	Left	238	201	80	73	208	P.54	XL	Left	216	189	79	76	215	P.64		
XL	Right	243	210	84	77	213	P.54	XL	Right	230	199	80	74	201	P.64		
XL1	Left	233	204	85	73	238	P.56	L	Left	204	178	70	67	208	P.82		
XL1	Right	230	195	78	77	235	P.56	L	Right	208	180	71	67	198	P.82		
L	Left	234	205	75	76	219	P.58	L1	Left	199	167	65	73	203	P.70		
L	Right	240	207	80	75	208	P.58	L1	Right	203	168	65	71	203	P.70		
L1	Left	214	192	72	66	211	P.60	L2	Left	204	180	73	70	213	P.84		
L1	Right	226	193	78	65	208	P.60	L2	Right	203	179	71	67	203	P.84		
L2	Left	228	185	71	76	240	P.62	LM	Left	189	167	63	66	200	P.76		
L2	Right	229	196	75	77	240	P.62	LM	Right	198	163	64	69	196	P.76		
LM	Left	216	189	79	76	215	P.64	LM1	Left	190	165	61	69	195	P.78		
LM	Right	230	199	80	74	201	P.64	LM1	Right	200	161	64	70	198	P.78		
M	Left	210	185	72	68	213	P.66	M	Left	189	174	68	69	203	P.86		
M	Right	215	188	71	69	200	P.66	M	Right	185	167	66	68	198	P.86		
M1	Left	207	185	69	67	198	P.68	MS	Left	190	169	64	70	193	P.80		
M1	Right	213	180	72	67	198	P.68	MS	Right	193	167	64	70	195	P.80		
M2	Left	199	167	65	73	203	P.70	MS2	Left	183	163	61	65	206	P.90		
M2	Right	203	168	65	71	203	P.70	MS2	Right	185	164	64	64	206	P.90		
M3	Left	216	173	67	72	213	P.72	S	Left	180	159	60	60	197	P.88		
M3	Right	204	166	69	67	223	P.72	S	Right	186	165	62	62	195	P.88		
M4	Left	210	179	71	70	210	P.74	S1	Left	175	145	55	65	179	P.92		
M4	Right	207	173	68	67	203	P.74	S1	Right	175	150	55	67	185	P.92		
MS	Left	189	167	63	66	200	P.76										
MS	Right	198	163	64	69	196	P.76										
MS1	Left	190	165	61	69	195	P.78										
MS1	Right	200	161	64	70	198	P.78										
S	Left	190	169	64	70	193	P.80										
S	Right	193	167	64	70	195	P.80										
CL	Left	190	169	64	70	193	P.80	CL	Left	175	145	55	65	179	P.92		
CL	Right	193	167	64	70	195	P.80	CL	Right	175	150	55	67	185	P.92		
CL1	Left	183	163	61	65	206	P.90										
CL1	Right	185	164	64	64	206	P.90										
CM	Left	175	145	55	65	179	P.92	CM	Left	173	150	56	63	175	P.94		
CM	Right	175	150	55	67	185	P.92	CM	Right	166	139	55	64	180	P.94		
CS	Left	140	119	51	48	165	P.96	CS	Left	140	119	51	48	165	P.96		
CS	Right	143	123	51	49	155	P.96	CS	Right	143	123	51	49	155	P.96		
CXS	Left	125	115	43	45	165	P.98	CXS	Left	125	115	43	45	165	P.98		
CXS	Right	130	119	45	43	158	P.98	CXS	Right	130	119	45	43	158	P.98		

The above information is for reference only. Please contact our technicians if you are not sure which size to choose.