

Silicone Prosthetic toe

Model TOE



Product descriptions:

- High definition silicone prosthetic toe(s)
- Design for toe amputation

Silicone TOE series

Number of Toes amputated	Recommended Models	Footband (Standard Footband is set at the metatarsal)
Single Toe	TOE1S	Without Footband
Single Toe	TOE1C, TOE2C, TOE3C, TOE4C, TOE5C	With Footband
Two Toes	TOE12C, TOE13C, TOE14C, TOE15C, TOE23C, TOE24C, TOE25C, TOE34C, TOE35C, TOE45C	With Footband
Three Toes	TOE123C, TOE124C, TOE125C, TOE134C, TOE135C, TOE145C, TOE234C, TOE235C, TOE245C, TOE345C	With Footband
Four or Five Toes	HDSFHT, HDSF	-

How to Order :

For Model TOE Series toe amputation: **Order form S5**
(All TOE series must be ordered in Semi-Custom Made (SCM) and with Silicone Filling)



TOE1S – Great Toe without Footband

TOE1C – Great Toe with Footband



TOE2C – Second Toe with Footband



TOE23C – Second and Third Toe with Footband



TOE34C - Third and Forth Toe with Footband
(Customer may request to cut the Footband to desired shape and length)



Without filling



Silicone filling